



Europäisches Patentamt
European Patent Office
Office européen des brevets



Publication number: **0 345 549 B1**

⑫

EUROPEAN PATENT SPECIFICATION

- ⑭ Date of publication of patent specification: 26.04.95 ⑮ Int. Cl.⁶: **C07K 1/18, C07K 16/00, G01N 33/68, //A61K39/395, C12P21/00**
- ⑰ Application number: 89109455.9
- ⑱ Date of filing: 25.05.89

⑤④ Removal of protein A from antibody preparations.

③① Priority: 08.06.88 US 204054

④③ Date of publication of application:
13.12.89 Bulletin 89/50

④⑤ Publication of the grant of the patent:
26.04.95 Bulletin 95/17

⑥④ Designated Contracting States:
AT BE CH DE ES FR GB GR IT LI LU NL SE

⑤⑥ References cited:

JOURNAL OF IMMUNOLOGICAL METHODS,
vol. 83, 1985, pages 169-177, Elsevier Science
Publishers B.V. (Biomedical Division); M.T.
DERTZBAUGH et al.: "An enzyme immunoassay
for the detection of staphylococcal protein
A in affinity-purified products"

⑦③ Proprietor: MILES INC.
1127 Myrtle Street
Elkhart
Indiana 46514 (US)

⑦② Inventor: Bloom, James W.
3663 May Road
Richmond
California, 94803 (US)
Inventor: Wong, Melvin F.
4324 25th Street
San Francisco
California 94114 (US)
Inventor: Mitra, Gautam
40 Cowper Avenue
Kensington
California, 94707 (US)

⑦④ Representative: Adrian, Albert, Dr. et al
Bayer AG
Konzernverwaltung RP
Lizenzen und Technische Kooperationen
D-51368 Leverkusen (DE)

Note: Within nine months from the publication of the mention of the grant of the European patent, any person may give notice to the European Patent Office of opposition to the European patent granted. Notice of opposition shall be filed in a written reasoned statement. It shall not be deemed to have been filed until the opposition fee has been paid (Art. 99(1) European patent convention).

CANCER RESEARCH, vol. 44, February 1984, pages 734-743; J. BALINT, Jr. et al.: "Tumoricidal response following perfusion over immobilized protein A: Identification of immunoglobulin oligomers in serum after perfusion and their partial characterization"

JOURNAL OF IMMUNOLOGICAL METHODS, vol. 93, 1986, pages 63-70, Elsevier Science Publishers B.V. (Biomedical Division); A. WARNES et al.: "Development of an enzyme-linked immunosorbent assay for staphylococcal protein A produced in Escherichia coli by pUC8 based plasmids containing the Staphylococcus aureus cowan I protein A gene"

METHODS OF BIOMEDICAL ANALYSIS, vol. 26, 1980, pages 1-45; E.E. BAYER et al.: "The use of the avidin-biotin complex as a tool in molecular biology"

JOURNAL OF IMMUNOLOGICAL METHODS, vol. 47, 1981, pages 129-144, Elsevier/North-Holland Biomedical Press; A.G. FARR et al.: "Immunohistochemistry with enzyme-labeled antibodies: a brief review"

FEBS LETTERS, vol. 28, no. 1, November 1972, pages 73-76, North-Holland Publishing Co., Amsterdam, NL; H. HJELM et al.: "Protein A from Staphylococcus aureus. Its isolation by affinity chromatography and its use as an immunosorbent for isolation of immunoglobulins"

JOURNAL OF IMMUNOLOGICAL METHODS, vol. 55, 1982, pages 277-296, Elsevier Biomedical Press; J.J. LANGONE: "Applications of immobilized protein A in

immunochemical techniques"

JOURNAL OF IMMUNOLOGICAL METHODS, vol. 62, 1983, pages 1-13, Elsevier Science Publishers B.V.; R. LINDMARK et al.: "Binding of immunoglobulins to protein A and immunoglobulin levels in mammalian sera"

Description

This disclosure is concerned generally with the purification of biologically active proteins and specifically with the reduction of Protein A contamination in therapeutic antibody preparations.

5 Protein A is a well known substance obtained from Staphylococcus aureus, Cowan Strain I. It has a high degree of antibody specificity and the protein has long been used to complex with antibodies. Protein A is commonly used in an immobilized form by being attached to various water insoluble support materials. The immobilized protein A is then used to complex with soluble antibodies which are subsequently eluted from antibody - immobilized protein A complexes. Protein A chromatography is thus an excellent method for
10 purifying ascites or tissue culture fluid derived monoclonal antibodies to homogeneity because of its simplicity and high degree of antibody specificity.

If the purified monoclonal antibodies are intended for therapeutic usage, however, a major safety concern is the possible presence of solubilized Protein A in the purified therapeutic product. Such solubilized Protein A is thought to result from the unintended detachment of Protein A from its support
15 material during the purification process.

Numerous publications link Protein A with toxicity and mitogenicity in animal models and humans (see, for example, Bensinger et al., J. Biol. Resp. Modif. 3, 347, 1984; Messerschmidt et al., J. Biol. Resp. Modif. 3,325, 1984; Terman and Bertram, Eur. J. Cancer Clin. Oncol. 21, 1115; 1985; and Ventura et al., Hortobagyl. Cancer Treat Rep. 71,411, 1987).

20 Unfortunately, to date there have been no assay methods available to measure very low amounts of Protein A (i.e., less than 15 pg of Protein A per mg of protein) that might be undesirably present in an antibody preparation intended for therapeutic use. Surprisingly, we found such an assay is now possible. Our assay led to a method of reducing Protein A contamination in antibody preparations to very low levels. The method thus permits the production of therapeutic antibody preparations using immobilized protein A
25 while assuring a low degree of Protein A contamination. Details of our assay and purification methods are described below.

To evaluate the amount of Protein A column leakage, if any, into monoclonal antibodies purified with immobilized Protein A we 1) developed a Protein A ELISA sensitive to the subnanogram range; 2) used the ELISA to determine Protein A levels in monoclonal antibodies purified by Protein A chromatography; and 3)
30 devised a method to reduce solubilized Protein A found to be contaminating the antibody preparations. Our highly sensitive Protein A ELISA was made possible by using biotinylated anti-Protein A as the immunoassay label, thus allowing measurement of less than 15 ng of Protein A per mg of antibody. Our method of reducing the amount of Protein A in an antibody - Protein A mixture comprises contacting the mixture with an anion exchange column under conditions sufficient to complex both components of the mixture onto an
35 anion exchange material and then selectively eluting the components by carefully varying ionic strength under conditions sufficient to assure elution of an antibody preparation substantially free of Protein A (i.e. less than 15 ng of Protein A per mg of antibody). Preferably the Protein A content is less than about 1 ng per mg of antibody, and the purification step is accomplished by applying the mixture to an ion exchange column such as a DEAE Trisacryl M or DEAE Sepharose column and then eluting the antibody with a NaCl
40 solution concentration gradient of about 0.025M to 0.25M as described below.

BRIEF DESCRIPTION OF THE FIGURES

Figure 1 is a graph showing an ELISA for Protein A at concentrations ranging from 0.02 mg/ml to
45 approximately 10 ng/ml.

Figure 2 is a graph showing a redrawn portion of the standard curve from 0.02 ng/ml to 1.3 mg/ml.

Figure 3 is a graph showing the column profile for the elution of both antibody and Protein A under conditions of increasing ionic strength.

50 SPECIFIC EMBODIMENTS

An enzyme labeled immunosorbent assay (ELISA) for low levels of Protein A was devised by biotinylating anti-Protein A IgG as follows:

Reagents: FTA Hemagglutination buffer (PBS) was obtained from BBL, Microbiology Systems (Cockeysville, MD). TMB (tetramethyl benzidine) Microwell Peroxidase Substrate and Peroxidase Substrate Solution B
55 were obtained from Kirkegaard & Perry Labs, Inc. (Gaithersburg, MD). Rabbit anti-Protein A IgG (#45F-8806) and bovine albumin were obtained from Sigma (St. Louis, MO). HRP-Streptavidin was purchased from Zymed Laboratories, Inc. (San Francisco, CA) and Biotin-X-NHS from Calbiochem (La Jolla, CA). Natural

Protein A AvidGel F was obtained from BioProbe International, Inc. (Tustin, CA). Recombinant Protein A was obtained from Repligen Corp. (Cambridge, MA). DEAE-Trisacryl M was obtained from LKB Instruments, Inc. (Gaithersburg, MD).

Methods: Rabbit anti-Protein A IgG was biotinylated according to the supplier's instructions. The procedure consisted of addition of 20 μ l of 20 mg/ml Biotin-X-NHS to a 1 mg/ml antibody solution in 0.1 M bicarbonate buffer followed by gentle agitation for 1 hour at room temperature. Excess Biotin-X-NHS was then removed by dialysis against PBS buffer.

The concentration of the Protein A standard was determined using an extinction coefficient of 1.46 for a 1% solution at 275 nm (see Sjoquist et al., Eur. J. Biochem 29,572, 1972)

Our ELISA (Enzyme labeled immunosorbant assay) was performed as follows: Nunc Immuno Plate (#4-3945, InterLab, Thousand Oaks, CA) 96 well flat bottom ELISA plates were coated with 100 μ l of 5 μ g/ml anti-Protein A IgG diluted into 0.05M carbonate, pH 9.6 and incubated overnight at 5 C. The plates were then washed with wash buffer (PBS-0.05% Tween 20), blocked with 200 μ l of wash buffer plus 1.5% BSA per well and incubated for one hour at 37 C. The plates were washed. Samples and the standard diluted in the wash buffer plus BSA were added in 100 μ l aliquots to the plate and incubated for one hour at 37 C. Wells were set aside that contained only 100 μ l of diluent to serve as blanks. The samples were discarded and the plates washed again with wash buffer. To each well 100 μ l of biotinylated anti-Protein A IgG diluted in wash buffer plus BSA (approximately 1/16,000) was applied and incubated for 1 hr at 37 C. The antibody was discarded and the plates washed. Finally, 100 μ l per well of Streptavidin-HRP was added and the plates were incubated for 1 hr at 37 C. The plates were washed and 100 μ l of a 1:1 solution of TMB Microwell Peroxidase Substrate and Peroxidase Substrate Solution B was added per well and incubated for approximately 5 minutes at room temperature. 100 μ l of 1N HCl was added to stop the reaction. The intensity of yellow color generated was proportional to the amount of Protein A present and was determined by measuring the absorbance at 450 nm on a Dynatech MR600 microtiter plate reader (Dynatech, Burlington, MA). The reading for each well was compensated for non-specific contributions by dividing by the absorbance at 570 nm. The average value of the blanks was subtracted to eliminate background. Experiments were done in duplicate on the same ELISA plate and average values utilized in the data analysis method.

ELISA RESULTS

ELISA sensitivity:

Purified Protein A was tested in the ELISA at concentrations ranging from 0.02 ng/ml to 10 ng/ml. The results shown in Figure 1 indicate saturation at approximately 10 ng/ml. The portion of the standard curve from 0.02 ng/ml to 1.3 ng/ml is redrawn in Figure 2 with expanded scales and shows that the linear portion of the curve is approximately in the range of 0.05 ng/ml to 0.6 ng/ml.

ELISA validation:

The Protein A ELISA was developed to quantitate Protein A in the presence of monoclonal antibodies. To determine whether the presence of a monoclonal antibody in stoichiometric excess would interfere with the assay the following experiment was undertaken. A 1 mg/ml solution of a monoclonal antibody purified to homogeneity by ion exchange chromatography was spiked with 10 μ g/ml of purified Protein A. As a control and standard the ELISA dilution buffer was also spiked with 10 μ g/ml Protein A. The mixtures were incubated for 1/2 hour at room temperature and then diluted to 20 ng/ml in the assay dilution buffer for assay. The results are shown in Figure 2. The two curves are nearly superimposable. The results for the standard and spiked monoclonal antibody were 20 ng/ml and 18.3 ng/ml respectively. The difference of 9% can be attributed to experimental error. These results suggest that the presence of an excess amount of a monoclonal antibody does not significantly interfere with the ELISA.

Leakage of Protein A:

The Protein A contamination levels of several monoclonal antibody preparations purified by Protein A chromatography were determined (Table I). Eluate pools 1a to 1d represent back-to-back 0.8 liter Protein A column runs with a 3 column volume purge of 0.2M glycine, pH 2.8 between runs. An initial Protein A level of approximately 300 ng/ml dropped to approximately 40 - 100 ng/ml in subsequent runs.

EP 0 345 549 B1

Eluate pools 2 and 3 represent a 1.5 liter Protein A column containing a blend of the immobilized Protein A used to obtain Pool 1 eluates and new immobilized Protein A. Pool 2 eluate Protein A values are similar to Pool 1 but Pool 3 eluate values are considerably less. The amount of Protein A leakage seems to decrease with column usage.

TABLE I

PROTEIN A LEAKAGE FROM IMMOBILIZED PROTEIN A COLUMNS USED TO PURIFY A MONOCLONAL ANTIBODY		
Eluate Pool #	Protein A(ng/ml)	Protein A/Monoclonal(ng/mg)
1a	316.8	630
1b	40.9	61
1c	87.8	114
1d	108.4	131
2a	316.0	238
2b	67.4	46
3a	19.6	28
3b	9.5	11

EXAMPLES

(Monoclonal Antibody Purification)

As measurable amounts of Protein A were present in monoclonal antibodies purified by immobilized Protein A, we looked for a method to reduce the Protein A contaminant levels. In an initial experiment, a monoclonal antibody previously purified by Protein A chromatography was spiked with purified Protein A to a concentration of 0.59 ug Protein A per mg monoclonal antibody and subjected to DEAE chromatography. The column profile is shown in Figure 3. Protein A eluted at a higher NaCl concentration than the IgG and good separation of the protein elution peaks was achieved.

To evaluate further the ability of DEAE chromatography to reduce Protein A contamination in a monoclonal antibody purified on a Protein A column, a Protein A eluate was chromatographed on DEAE (see Table 2).

TABLE 2

DEAE CHROMATOGRAPHY OF A MONOCLONAL ANTIBODY PURIFIED ON IMMOBILIZED PROTEIN A			
Fraction	IgG(mg/ml)	Protein A(ng/ml)	Protein A/IgG(ng/mg)
Leading edge	0.69	1.0	1.45
Main peak	3.13	7.8	2.49
Trailing edge	0.45	87.6	194.67

Significant separation of Protein A from the monoclonal was observed in confirmation of the spiking experiment. Details of a representative purification process are shown in the following example.

Monoclonal antibody purification: a IgG1 tissue culture fluid derived anti coagulant Factor VIII monoclonal antibody (designated C7F7) was purified in the following manner: 1) Tissue culture fluid is clarified by filtration, 2) Polyethylene glycol is added (17% w/v) to the supernatant solution and dissolved. Precipitate is separated and the solution is discarded. The precipitate is dissolved in 0.05 M tris-(hydroxymethyl)aminomethane, 0.15 M sodium chloride, pH 8.00, to a volume of 2.5% of the original tissue culture fluid. Dissolved precipitate may be stored frozen at -20 C or colder; 3) Fresh or thawed dissolved precipitate is clarified by centrifugation and/or filtration; 4) Solution is contacted with Protein A Avidgel F (R) (or equivalent immobilized Protein A) equilibrated with dissolving buffer and washed with same. C7F7 is removed by elution with 0.05 M sodium acetate, 0.15 M sodium chloride pH 4.00; 5). The eluate is diafiltered against not less than 6 volumes of DEAE equilibration buffer (0.025 M sodium chloride, 0.025 M

tris(hydroxymethyl)aminomethane, pH 8.60); 6) The solution is contacted with DEAE Sepharose (R) or equivalent anion exchange resin (previously equilibrated) and washed with DEAE equilibration buffer. C7F7 is eluted with a sodium chloride gradient from 0.025 M to 0.25 M. Eluate is collected based on A280. The trailing edge is discarded. The trailing edge is defined as less than 20% of the maximum peak A280; 7) Diafilter solution against not less than 6 volumes of phosphate buffered saline (8.9 mM disodium phosphate, 0.7 mM monosodium phosphate, 1.6 mM monopotassium phosphate, 0.15 M sodium chloride). Using the above methods, the Protein A content of the antibody preparation was reduced to less than 15 ng/mg antibody (i.e. range of Protein A was from 0.9 to 14 ng/mg of antibody).

Protein A levels in several Protein A purified monoclonal antibody preparations further purified with DEAE chromatography were found to be in the range 0.9 to 14 ng/mg.

DISCUSSION

Numerous ELISAs have been developed for the detection of Protein A in the nanogram range of sensitivity (Maxim et al., J. Clin. Microbiol. 4,418, 1976; Langone et al., J. Immunol. Methods 18,281, 1977; Fey and Burkhard, J. Immunol. Methods 47, 99, 1981; Lofdahl et al., Proc. Nat. Acad. Sci., U.S.A., 80,697, 1983; Olsvik and Berdal, Acta Pathol. Microbiol. Immunol. Scand. Sect. B. Microbiol. 89, 289, 1981; Dertzbaugh et al., J. Immunol. methods 83, 169, 1985; Considine et al., Bios. Rep. 6, 933, 1986; Warnes, et al., J. Immunol. Methods, 93, 63, 1986). As our application involves assaying Protein A in the presence of excess amounts of IgG, the assay had to be capable of detecting both free Protein A and Protein A complexed to the Fc region of IgG. Thus the antibodies used in the ELISA have to be specific for epitopes on the Protein A molecule. The majority of the published Protein A ELISA techniques are unsuitable for this application because they utilize the Fc binding ability of Protein A in the assay (Maxim, et al., 1976; Langone et al., 1977; Fey and Burkhard, 1981; Lofdahl, et al., 1983; Olsvik and Berdal, 1981; Considine, et al., 1986; Warnes, et al., 1986). We have developed an ELISA utilizing an anti-Protein A coating antibody and biotinylated anti-Protein A as the detection antibody. The biotinylated secondary antibody ELISA had an assay sensitivity approximately one hundred fold greater than a similar assay using an alkaline phosphatase labeled reagent (Dertzbaugh et al., 1985). In addition, with a sensitivity less than 0.1 ng/ml, the ELISA is five to ten fold more sensitive than the most sensitive Protein A system yet published (Warnes, et al., 1986). Replacement of the rabbit anti-Protein A-biotin with rabbit anti-Protein A followed by goat anti-rabbit IgG-biotin might improve the assay sensitivity even more (Warnes, et al., 1986).

Protein A contamination of an affinity column purified monoclonal antibody has been previously reported (Dertzbaugh et al., J. Immunol. Methods 83, 169, 1985) and the Protein A values appear to be similar to those reported here. The 2-fluoro-1-methyl pyridinium toluene-4-sulfonate (FMP) activated gel used to link Protein A in our studies is reported to be a more stable linkage than the cyanogen bromide-activated gel (Ngo, Bio/Technology 4, 134, 1986) used by Dertzbaugh, et al. The observations that the leakage rates of the two Protein A matrices are similar and that usage seems to reduce leakage, support the hypothesis that the majority of the Protein A that leached off the columns was noncovalently bond to the matrix (Dertzbaugh et al., J. Immunol. Methods, 1985).

DEAE chromatography was demonstrated to effectively reduce Protein A contaminant levels in preparations of a mouse monoclonal antibody. In the presence of excess rabbit or human IgG, Protein A has been shown to form complexes with the molecular formula of [IgG2ProteinA] or the dimeric structure [(IgG)-2ProteinA]2 (Balint et al., Cancer Res. 44, 734, 1984; Das et al., Anal. Biochem., 145, 27, 1985). These complexes could be separated from IgG by Sepharose CL-6B (Balint et al., Cancer Res. 44, 734, 1984) or by gel filtration HPLC. Model complexes formed by adding protein A to serum or to monomeric IgG have been reported to activate Fc-bearing leukocytes and the complement system to generate oxidant and anaphylatoxin activity in vitro (Balint et al., Cancer Res. 44, 734, 1984). A purification step, such as DEAE or gel filtration chromatography, designed to reduce contaminant Protein A levels would minimize the possibility of toxic effects in recipients of therapeutic monoclonal antibodies.

Given the above disclosure and examples, it is thought that numerous variations will occur to those skilled in the art. Accordingly, it is intended that the above examples should be construed as illustrative and that the scope of the invention disclosed herein should be limited only by the following claims.

Claims

1. A method of reducing the amount of Protein A from a mixture of antibodies and Protein A to less than about 15 ng of Protein A per mg of antibodies, the method comprising the steps of contacting the mixture with an anion exchange material under conditions sufficient to adsorb both the antibodies and

the Protein A and then sequentially eluting the antibodies and Protein A under conditions of increasing ionic strength.

2. The method of claim 1 wherein the conditions of increasing ionic strength comprise eluting the antibody and Protein A with a NaCl solution having a concentration gradient of about 0.025 M to 0.25 M and the amount of Protein A is reduced to less than about 1 ng per mg of antibodies.
3. The method of claim 1 wherein the anion exchange material comprises Diethylaminoethyl (DEAE) Trisacryl M or DEAE Sepharose.
4. The method of claim 1 wherein the antibody preparation obtained is substantially free of Protein A.
5. The method of claim 1 wherein the antibody preparation comprises at least one type of monoclonal antibody.
6. The method of claim 1 wherein the antibody preparation comprises monoclonal antibodies to blood coagulation Factor VIII and the Protein A content is less than about 15 ng per mg of antibodies.
7. The method of claim 6 wherein the conditions of increasing ionic strength comprise eluting the antibody and Protein A with a NaCl solution having a concentration gradient of about 0.025 M to 0.25 M and the amount of Protein A is reduced to less than about 1 ng per mg of antibody.
8. An ELISA for determining the concentration of Protein A in an antibody solution comprising the steps of contacting the solution with a immobilized antibody that binds to the Protein A and using that binding in an immunoassay with labelled antibodies that bind to Protein A to determine the concentration of Protein A, the improvement comprising using biotinylated anti-Protein A antibodies as the labelled antibodies.

Patentansprüche

1. Verfahren zur Verringerung der Menge an Protein A aus einer Mischung von Antikörpern und Protein A auf weniger als ca. 15 ng Protein A pro mg Antikörper, wobei das Verfahren Stufen umfaßt, in denen man die Mischung mit einem Anionenaustauschmaterial unter Bedingungen zusammenbringt, die ausreichen, sowohl die Antikörper als auch das Protein A zu adsorbieren, und man dann die Antikörper und Protein A unter Bedingungen steigender Ionenstärke sequentiell eluiert.
2. Verfahren gemäß Anspruch 1, worin die Bedingungen steigender Ionenstärke darauf beruhen, daß der Antikörper und Protein A mit einer NaCl-Lösung, die einen Konzentrationsgradient von ca. 0,025 M bis 0,25 M aufweist, eluiert und die Menge von Protein A auf weniger als ca. 1 ng pro mg Antikörper verringert werden.
3. Verfahren gemäß Anspruch 1, worin das Anionenaustauschmaterial Diethylaminoethyl(DEAE)-Trisacryl M oder DEAE-Sepharose enthält.
4. Verfahren gemäß Anspruch 1, worin die erhaltene Antikörper-Zubereitung im wesentlichen frei von Protein A ist.
5. Verfahren gemäß Anspruch 1, worin die Antikörper-Zubereitung mindestens einen Typ eines monoklonalen Antikörpers enthält.
6. Verfahren gemäß Anspruch 1, worin die Antikörper-Zubereitung monoklonale Antikörper zu Blutkoagulationsfaktor VIII enthält und der Protein A-Gehalt weniger als ca. 15 ng pro mg Antikörper beträgt.
7. Verfahren gemäß Anspruch 6, worin die Bedingungen steigender Ionenstärke darauf beruhen, daß man den Antikörper und Protein A mit einer NaCl-Lösung, die einen Konzentrationsgradient von ca. 0,025 M bis 0,25 M aufweist, eluiert und die Menge von Protein A auf weniger als ca. 1 ng/mg Antikörper verringert ist.

8. ELISA zur Bestimmung der Konzentration von Protein A in einer Antikörper-Lösung, in welchem die Stufen enthalten sind, in denen man die Lösung mit einem immobilisierten Antikörper, der an das Protein A gebunden wird, zusammenbringt und man diese Bindung in einem Immunoassay mit markierten Antikörpern, die sich an Protein A binden, nutzt, um die Konzentration von Protein A zu bestimmen, wobei sich eine Verbesserung des Verfahrens dadurch ergibt, daß man biotinylierten Anti-Protein-A-Antikörper als den markierten Antikörper verwendet.

Revendications

1. Procédé pour réduire la quantité de Protéine A d'un mélange d'anticorps et de Protéine A à moins d'environ 15 ng de Protéine A par mg d'anticorps, ledit procédé comprenant les étapes de mise en contact du mélange avec une matière échangeuse d'anions dans des conditions suffisantes pour adsorber à la fois les anticorps et la Protéine A, puis d'élution séquentielle des anticorps et de la Protéine A dans des conditions de force ionique croissante.
2. Procédé suivant la revendication 1, dans lequel les conditions de force ionique croissante comprennent l'élution des anticorps et de la Protéine A avec une solution de NaCl ayant un gradient de concentration d'environ 0,025 M à 0,25 M, la quantité de Protéine A étant réduite à moins d'environ 1 ng/mg d'anticorps.
3. Procédé suivant la revendication 1, dans lequel la matière échangeuse d'anions consiste en diéthylamino-éthyl (DEAE)-Trisacryl M ou DEAE-Sepharose.
4. Procédé suivant la revendication 1, dans lequel la préparation d'anticorps obtenue est pratiquement dépourvue de Protéine A.
5. Procédé suivant la revendication 1, dans lequel la préparation d'anticorps comprend au moins un type d'anticorps monoclonal.
6. Procédé suivant la revendication 1, dans lequel la préparation d'anticorps comprend des anticorps monoclonaux contre le Facteur VIII de coagulation sanguine et la teneur en Protéine A est inférieure à environ 15 ng par mg d'anticorps.
7. Procédé suivant la revendication 6, dans lequel les conditions de force ionique croissante comprennent l'élution des anticorps et de la Protéine A avec une solution de NaCl ayant un gradient de concentration d'environ 0,025 M à 0,25 M, la quantité de Protéine A étant réduite à moins d'environ 1 ng par mg d'anticorps.
8. Analyse par immunosorbant lié à un enzyme (ELISA) pour la détermination de la concentration de Protéine A dans une solution d'anticorps, comprenant les étapes de mise en contact de la solution avec un anticorps immobilisé qui se lie à la Protéine A et l'utilisation de cette liaison dans une analyse immunologique avec des anticorps marqués qui se lient à la Protéine A pour déterminer la concentration de Protéine A, le perfectionnement comprenant l'utilisation d'anticorps anti-Protéine A biotinylés comme anticorps marqués.



